

Molecular Design and Synthesis

- > Links to Session 1: Light/Matter
- > History of Making Things
- > Structure and Reactivity 101
- > Design Process/Matter Lifecycle
- > Polymeric Metal Complexes

> Design Process/Matter Lifecycle

- Target molecule
- Synthetic strategy
- Chemicals
- Laboratory notebook
- Reaction setup
- Isolation
- Purification
- Then where?


> Design Process/Matter Lifecycle

- Target molecule
- Synthetic strategy
- Chemicals
- Laboratory notebook
- Reaction setup
- Isolation
- Purification
- Then where?


Traditional Synthesis

Combinatorial Synthesis


•Target molecule



Callipeltoside A



Callipeltoside B




Callipeltoside C

Isolation: From marine sponges, Minale et al, University of Naples.
2.5 kg of sponge furnished 3.5 mg of the natural product.

Activity: Callipeltoside A inhibit in vitro proliferation of cancer cells; HIV infected cell protection.

Previous Total Synthesis: Paterson, Evans, and Trost. **Fragment Synthesis:** Hoye, Olivo

From: Patrick Eidem, Marshall group, UVA Chemistry Department




From: Patrick Eidem, Marshall group, UVA Chemistry Department

• chemicals: selection




chemicals: selection



Tetrabromobisphenol A, see 33,039-6, 4,4'-Isopropylidenebis(2,6-dibromophenol) page 1102			
41,122-1 ★	2,2',6,6'-Tetrabromobisphenol A diallyl ether, 99% [25327-89-3] (CH ₂) ₂ C(C ₆ H ₂ (Br) ₂ OCH ₂ CH=CH ₂) ₂ FW 624.03 mp 118-120°	250g 1kg	34.50 97.80
43,747-6 ★	2,2',6,6'-Tetrabromobisphenol A ethoxylate (1 EO/phenol) diacrylate [66710-97-2] [H ₂ C=CHCO ₂ (CH ₂ CH ₂ O) _n C ₆ H ₂ (Br) ₂ -4- ₂ C(CH ₃) ₂ mp 136-138° IRRITANT Average M _n ca. 740. Contains 300-500 ppm monomethyl ether hydroquinone	100g 500g	26.00 86.50
7480-3	Tetrabromocatechol, 96% [488-47-1] (tetrabromopyrocatechol) C ₆ Br ₄ -1,2-(OH) ₂ FW 425.72 mp 189-193° Beil. 6,786 FT-NMR 1(2),298C FT-IR 1(1),1103C Safety 2.3232B R&S 1(1),1295A RTECS# UX2430000 IRRITANT	5g 25g	20.70 68.60
27,708-8 ★	3,4,5,6-Tetrabromo- <i>o</i> -cresol, 98% [576-55-6] CH ₃ C ₆ Br ₄ OH FW 423.75 mp 209-212° Beil. 6,362 Merck Index 13,9260 FT-NMR 1(2),293C Safety 2.3232C R&S 1(1),1291B RTECS# GP3135000 IRRITANT	100g 500g	34.30 111.40
3',3'',5',5''-Tetrabromo- <i>m</i> -cresolsulfonephthalein, see Bromocresol Green			

<chem>O=C(O)CN1CCN(CC1)CC(=O)O</chem> • 4HCl • xH ₂ O	<chem>O=C(O)c1ccccc1OC2C(=O)c3ccccc3OC2N=C=S</chem>	<chem>Oc1ccccc1COc2ccccc2OC3C(=O)c4ccccc4OC3O</chem>
40,040-8	33,562-2	44,741-2
<chem>Oc1ccccc1COc2ccccc2OC3C(=O)c4ccccc4OC3O</chem>	<chem>Oc1ccccc1COc2ccccc2OC3C(=O)c4ccccc4OC3O</chem>	<chem>Oc1ccccc1COc2ccccc2OC3C(=O)c4ccccc4OC3O</chem>

chemicals: ordering



• chemicals: organics




• chemicals: Inorganics



• chemicals: solvents



solvent cabinet



solvent stills



solvent columns

•Workspace: Bench




- workspace: Fume hood



- workspace: Glovebox



• Laboratory notebook



5/7

$$\text{Fe} \left(\text{C}_6\text{H}_4\text{Cl} \right)_3 \xrightarrow[\text{3. HN(C}_2\text{H}_5)_2]{\text{1. ODCB, 2. C}_6\text{H}_4\text{Et}} \text{Fe} \left(\text{C}_6\text{H}_4\text{Et} \right)_3$$

PREPARATION OF A DIBLOCK CO-POLYMER:
POLY-2-ETHYL-2-OXAZOLINE-POLY-2-UNDECYL-2-OXAZOLINE

STEP 1: SYNTHESIS OF THE FIRST BLOCK

$$\text{Fe} \left(\text{C}_6\text{H}_4\text{Cl} \right)_3 \xrightarrow[\text{C}_6\text{H}_4\text{Et}]{\text{ODCB}} \text{Fe} \left(\text{C}_6\text{H}_4\text{Et} \right)_3$$

COMPOUND	MW	P	AMT	MOL	EQUIV
Fe ^{II} INITIATOR	1105.16		12.8mg	1.16×10^{-5}	1
ODCB		1.36	.91ml		
2-ETHYL-2-OXAZOLINE	99.13	.982	.207g .210ml	2.08×10^{-3}	180 30/site

TO THE Fe^{II} INITIATOR WAS ADDED ODCB THEN 2-ETHYL-2-OXAZOLINE THE RN WAS A DARK PURPLE AND STIRRED AT 83°C FROM 3³⁰ PM - 9¹⁰ PM (6h). AT 9:30 PM THE INITIATOR DID NOT DISSOLVE! PUT TEMPERATURE UP TO 105°C AT C

(1st BLOCK PROPAGATED FOR 13h)

5/8 AT 10:30 AM SOLN IS RED-PURPLE & INITIATOR IS REQUIRED. NEED HIGH TEMPERATURES. A SMALL ALIQUOT WAS TAKEN OUT (~45%) AND WAS TERMINATED W DIPROPYL AMINE. ~~THE REST WAS TREA~~

• Reaction setup



• Isolation





Separatory funnel





Rotary evaporator


• Purification



• Characterization: ^1H NMR spectroscopy



• Characterization: UV/vis spectroscopy





Iron Tris(bpy) PLA

• characterization: thermal methods




thermal analysis: decomposition, melting/glass transition temperature
(also: elemental analysis)


• synthetic scheme: C₃-C₉ Fragment




PMB-DMA, CSA,
toluene, reflux, 3 h
79% two steps




Red-Al, THF,
0 °C, 24 h, 87%




D-(-)-DIPT, TIP,
TBHP, 4A MS,
CH₂Cl₂, -20 °C, 83%



1. Red-Al, THF,
-20 °C, 18 h
2. TBSOTf,
2,6-Lutidine, rt, 3 h
74% two steps




DIBAL-H, CH₂Cl₂
-78 °C to -50 °C, 30 h,
74%




Swern,

-78 °C, 3 h, 100%



From: Patrick Eidem, Marshall group, UVA Chemistry Department

• Reaction optimization: selectivity



Metal	P	Ligand	Yield	R : S
Li	TBDPS	None	63	1 : 1
MgBr	TBDPS	None	63	3 : 1
ZnBr	TBDPS	None	68	4 : 1
ZnBr	TBDPS	(1R, 2S) NME	73	13 : 1

From: Patrick Eidem, Marshall group, UVA Chemistry Department

-Then where? In the literature



*Graphic from American Chemical Society Website

-Then where? on the shelf



-Then where? chemical waste




-Then where? fumes



-Then where? Practical uses

- drug testing
- property screening
- many other applications


• New Approach: Combinatorial synthesis




• New Approach: Combinatorial synthesis





• Manual vs robots: costs involved





OR



• 1st Combinatorial chemical library



• High throughput screening



•Combinatorial synthesis Laboratory



• Robot and multi-well trays

